

Frequency of Tests

drvalbers.com

Daily (Days X-rays Taken)

Sensitometry and densitometry

- Keep films for 60 days
- Keep QA logs until next State inspection

Processor temperature check

If you process less than 10 films per week, do QA first day films are processed

For days that films are not processed, chart this on the QA log

Quarterly (Every 3 months)

Repeat Analysis

First 1000 films

- Keep paperwork until next State inspection

Date	Date Completed	Initials

Semiannually (Every 6 months)

Darkroom fog test

- Keep films until next State inspection

Date	Date Completed	Initials

Annually (Every 12 months)

Radiation Program Audit

- Keep paperwork until next State inspection

Date	Date Completed	Initials

Biennially (Every 24 months)

Screen-film contact

- Keep films until next State inspection

Screen-film-cassette speed match

- Keep films until next State inspection

Integrity of aprons, gloves and gonad shielding

- Keep films until next state inspection

X-ray machine performance evaluation

- Keep paperwork from x-ray company until next State inspection

Date	Date Completed	Initials

Sensitometer and densitometer calibration

- Keep paperwork from Acurad Technical Services until next State inspection

Date	Date Completed	Initials

QA Records

- Keep all paperwork until next state inspection

- Personnel monitoring records must be preserved for the lifetime of the individual or a minimum of 30 years after termination of employment with a facility, whichever is less.